

## Product datasheet for **TP302833L**

### SIRT6 (NM\_016539) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human sirtuin (silent mating type information regulation 2 homolog) 6 ( <i>S. cerevisiae</i> ) (SIRT6), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202833 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MSVNYAAGLSPYADKGGKGLPEIFDPPEELERKVVWELARLVWQSSSVFHTGAGISTASGIPDFRGPHGV  
WTMEERGLAPKFDTTFESARPTQTHMALVQLERVGLLRFLVSNVDGLHVRSGFPRDKLAELHGNMFVEE  
CAKCKTQYVRDVTVMGLKATGRLCTVAKARGLRACRGELRDTILDWEDSLPDRDLALADEASRNADLS  
ITLGTSLQIRPSGNLPLATKRRGGRLVIVNLQPTKHDRHADLRHGYVDEVMTRLMKHLGLEIPAWDGP  
VLERALPPLPRPPTPKLEPKESPTRINGSIPAGPKQEPCAQHNGSEPAASPKRERPTSPAPHRPPKRVKA  
KAVPS

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	38.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_057623</a>




[View online »](#)

Locus ID: 51548  
UniProt ID: [Q8N6T7](#)  
RefSeq Size: 1657  
Cytogenetics: 19p13.3  
RefSeq ORF: 1065  
Synonyms: SIR2L6

**Summary:** This gene encodes a member of the sirtuin family of NAD-dependent enzymes that are implicated in cellular stress resistance, genomic stability, aging and energy homeostasis. The encoded protein is localized to the nucleus, exhibits ADP-ribosyl transferase and histone deacetylase activities, and plays a role in DNA repair, maintenance of telomeric chromatin, inflammation, lipid and glucose metabolism. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Mar 2016]

**Protein Families:** Druggable Genome, Transcription Factors

### Product images:



Coomassie blue staining of purified SIRT6 protein (Cat# [TP302833]). The protein was produced from HEK293T cells transfected with SIRT6 cDNA clone (Cat# [RC202833]) using MegaTran 2.0 (Cat# [TT210002]).