

Product datasheet for **TP302925L**

PLTP (NM_182676) Human Recombinant Protein

Product data:

| | |
|---------------------------------------|---|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant protein of human phospholipid transfer protein (PLTP), transcript variant 2, 1 mg |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | >RC202925 protein sequence Red =Cloning site Green =Tags(s) |

MALFGALFLALLAGAHAEFPGCKIRVTSKALELVKQEGLRFLEQELETITIPDLRGKEGHFYNNISEVKV
TELQLTSSSELDQPPQELMLQITNASLGLRFRRLLYWFLKVDLSTFITSGMRFLNQQICPVLYHAG
TVLLNSLLDTPVRSSVDELVGIDYSLMKDPVASTSNLDMDFRGAFFPLTERNWSPNRAVEPQLQEEER
MVYVAFSEFFFDSSAMESYFRAGALQLLLVGDKVPDMLLRATYFGSIVLLSPAVIDSPLKLELRVLAP
PRCTIKPSGTTISVTASVTIALVPPDQPEVQLSSMTMDARLSAKMALRGKALRTQLDLRRFRIYSNHSAL
ESLALIPLQAPLKTMLQIGVMPMLNERTWRGVQIPLPEGINFVHEVVTNHAGFLTIGADLHFAKGLREVI
EKNRPADVRASTAPTSTAAV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

| | |
|----------------|--|
| Tag: | C-Myc/DDK |
| Predicted MW: | 47.6 kDa |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol |
| Preparation: | Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps. |
| Note: | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage: | Store at -80°C. |
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |



[View online »](#)

RefSeq: [NP_872617](#)

Locus ID: 5360

UniProt ID: [P55058](#), [P55058-2](#)

RefSeq Size: 1944

Cytogenetics: 20q13.12

RefSeq ORF: 1323


Synonyms: BPIFE; HDLCQ9

Summary: The protein encoded by this gene is one of at least two lipid transfer proteins found in human plasma. The encoded protein transfers phospholipids from triglyceride-rich lipoproteins to high density lipoprotein (HDL). In addition to regulating the size of HDL particles, this protein may be involved in cholesterol metabolism. At least two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Secreted Protein

Protein Pathways: PPAR signaling pathway

Product images:



Coomassie blue staining of purified PLTP protein (Cat# [TP302925]). The protein was produced from HEK293T cells transfected with PLTP cDNA clone (Cat# [RC202925]) using MegaTran 2.0 (Cat# [TT210002]).