

## Product datasheet for **TP302937L**

### NRIP3 (NM\_020645) Human Recombinant Protein

#### Product data:

**Product Type:** Recombinant Proteins

**Description:** Recombinant protein of human nuclear receptor interacting protein 3 (NRIP3), 1 mg

**Species:** Human

**Expression Host:** HEK293T

**Expression cDNA Clone  
or AA Sequence:** >RC202937 protein sequence  
**Red**=Cloning site **Green**=Tags(s)

MFYSGLLTEGGRKETDMREAASLRQQRRMKQAVQFIHKDSADLLPLDGLKKGSSKDMQPHNILQRRRLME  
TNLSKLRSGPRVPWASKTNKLNQAKSEGLKKSEEDDMILVSCQCAGKDVKALVDTGCLYNLISLACVDRL  
GLKEHVKSHKHEGEKLSLPRHLKVVGGQIEHLVITLGSRLDCPAAWDDNEKNLSLGLQTLRSLKCIINL  
DKHRLIMGKTDKEEIPFVETVSLNEDNTSEA

**TR**TRPLE**QKLISEEDLAANDILDYKDDDDKV**

**Tag:** C-Myc/DDK

**Predicted MW:** 26.8 kDa

**Concentration:** >0.05 µg/µL as determined by microplate BCA method

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining

**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

**Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

**Note:** For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

**Storage:** Store at -80°C.

**Stability:** Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

**RefSeq:** [NP\\_065696](#)


**Locus ID:** 56675

**UniProt ID:** [Q9NQ35](#)



[View online »](#)

RefSeq Size:	3809
Cytogenetics:	11p15.4
RefSeq ORF:	723
Synonyms:	C11orf14; NY-SAR-105
Protein Families:	Druggable Genome

**Product images:**

Coomassie blue staining of purified NRIP3 protein (Cat# [TP302937]). The protein was produced from HEK293T cells transfected with NRIP3 cDNA clone (Cat# [RC202937]) using MegaTran 2.0 (Cat# [TT210002]).