

Product datasheet for **TP302969L**

PGAP4 (NM_032342) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 9 open reading frame 125 (C9orf125), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202969 protein sequence Red =Cloning site Green =Tags(s) MSTSTSPAAMLLRRLRRLSWGSTAVQLFILT VVTFGLLAPLACHRLLHSYFYLRHWHNLNQMSQEFLQQSL KEGEALHYFEELPSANGSVPIVWQATPRPWLVIITVDRQPGFHYVLQVVSQFHRLLQQCGPQCEGHQ LFLCNVERSVSHFDAKLLSKYVPVANRYEGTEDDYGDDPSTNSFEKEKQDYVYCLESSLQTYNPDYVLMV EDDAVPEEQIFPVLEHLLRARFSEPHLRDALYKLYHPERLQHYINPEPMRILEWVGVMLLGPLLTWIIY MRFASRPGFSWPVMLFFSLYSMGLVELVGRHYFLELRLSPSLYSVVPASQCCTPAMLFAPAARRTLTY LSQVYCHKGFGKDMALYSLLRAKGERAYVVEPNLVKHIGLFSSLRYNFHPSSL TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	46.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_115718</u>
Locus ID:	84302


[View online »](#)

UniProt ID:	<u>Q9BRR3</u>
RefSeq Size:	2110
Cytogenetics:	9q31.1
RefSeq ORF:	1209
Synonyms:	C9orf125; TMEM246
Summary:	Golgi-resident glycosylphosphatidylinositol (GPI)-N-acetylgalactosamine transferase involved in the lipid remodeling steps of GPI-anchor maturation. Lipid remodeling steps consist in the generation of 2 saturated fatty chains at the sn-2 position of GPI-anchors proteins (PubMed:29374258). Required for the initial step of GPI-GalNAc biosynthesis, transfers GalNAc to GPI in the Golgi after fatty acid remodeling by PGAP2 (PubMed:29374258). [UniProtKB/Swiss-Prot Function]
Protein Families:	Transmembrane

Product images:



Coomassie blue staining of purified TMEM246 protein (Cat# [TP302969]). The protein was produced from HEK293T cells transfected with TMEM246 cDNA clone (Cat# [RC202969]) using MegaTran 2.0 (Cat# [TT210002]).