

Product datasheet for **TP302972**

PRADC1 (NM_032319) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 2 open reading frame 7 (C2orf7), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202972 protein sequence Red =Cloning site Green =Tags(s)
	 MVPGAAGWCCLVLWLPACVAAHGFRIHDYLYFQVLSPGDIRYIFTATPAKDFGGIFHTRYEQIHLVPAEP PEACGELSNFFIQDQIALVERGGCSFLSKTRVVQEHGGRAVIISDनावDNDNSFYVEMIQDSTQRTADIP ALFLLGRDGYMIRRSLEQHGLPWAIISIPVNVTSIPTFELLQPPWTFW TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	20.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_115695</u>
Locus ID:	84279
UniProt ID:	<u>Q9BSG0</u>
RefSeq Size:	1105



[View online »](#)

Cytogenetics:	2p13.2
RefSeq ORF:	564
Synonyms:	C2orf7; PAP21
Summary:	Plays a role in the modulation of physical activity and adiposity.[UniProtKB/Swiss-Prot Function]
Protein Families:	Protease, Transmembrane

Product images:

Coomassie blue staining of purified PRADC1 protein (Cat# TP302972). The protein was produced from HEK293T cells transfected with PRADC1 cDNA clone (Cat# [RC202972]) using MegaTran 2.0 (Cat# [TT210002]).