

Product datasheet for **TP302976**

CDKN2AIPNL (NM_080656) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human CDKN2A interacting protein N-terminal like (CDKN2AIPNL), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC202976 protein sequence Red =Cloning site Green =Tags(s)
	MVGGEAAAVEELVSGVRQAADF AEQFRSYSESEKQWKARMEFILRHLPDYRDPPDGSGRLDQLLSLSMV WANHLFLGCSYNKDLLDKVMEMADGIEVEDLPQFTTRSELMKKHQS
	TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	13 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_542387</u>
Locus ID:	91368
UniProt ID:	<u>Q96HQ2</u>
RefSeq Size:	1270




[View online »](#)

Cytogenetics: 5q31.1

RefSeq ORF: 348

Synonyms: C2AIL

Product images:


Coomassie blue staining of purified CDKN2AIPNL protein (Cat# TP302976). The protein was produced from HEK293T cells transfected with CDKN2AIPNL cDNA clone (Cat# [RC202976]) using MegaTran 2.0 (Cat# [TT210002]).