

Product datasheet for TP303012L

TMEM128 (NM_032927) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human transmembrane protein 128 (TMEM128), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203012 protein sequence Red =Cloning site Green =Tags(s) MPPFGYKTETSTAVEKKEKPLPRLNIHSGFWILASIVVTTYVDFFKTLKENFHTSSWFLCGSALLLVSL IAFYCIVYLEWYCGIGEYDVKYPALIPITTASFIAAGICFNIALWHVWSFFTPLLLFTQFMGVVMFITLL G TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	15.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_116316
Locus ID:	85013
UniProt ID:	Q5BJH2
RefSeq Size:	1785


[View online »](#)

Cytogenetics: 4p16.3

RefSeq ORF: 423

Protein Families: Transmembrane

Product images:


Coomassie blue staining of purified TMEM128 protein (Cat# [TP303012]). The protein was produced from HEK293T cells transfected with TMEM128 cDNA clone (Cat# [RC203012]) using MegaTran 2.0 (Cat# [TT210002]).