

Product datasheet for TP303030M

Apc5 (ANAPC5) (NM_016237) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human anaphase promoting complex subunit 5 (ANAPC5), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203030 protein sequence Red=Cloning site Green=Tags(s)

MASVHESLYFNPMMTNGVVHANVFGIKDWVTPYKIAVLVLLNEMSRGTGEGAVSLMERRRLNQLLLPLLQG
PDITLSKLYKLIEESCPQLANSVQIRIKLMAEGELKDMEQFFDDLSDSFSGTEPEVHKTSVVGLFLRHMI
LAYSKLSFSQVFKLYTALQQYFQNGEKKTVEDADMELTSRDEGERKMEKEELDVSVREEEVSCSGPLSQK
QAEFFLSQQASLLKNDETKALTPASLQKELNNLLKFNPDFAEAHYLSYLNLRVQDVFSTHSLHLYFDR
LILTGAESKSNGEEGYGRSLRYAALNLAALHCRFGHYQQAELALQEAIRAQESNDHVCLQHCLSWLYVL
GQKRSDSYVLLHESVKKAVHFGLPYLASLGIQSLVQQRAFAGKTANKLMDALKDSDLHWHKHSSELIDI
SIAQKTAIWRLYGRSTMALQQAQMLLSMNSLEAVNAGVQQNNTESFAVALCHLAELHAEQGCFAAAASEVL
KHLKERFPNSQHAQLWMLCDQKIQFDRAMNDGKYHLADSLVTGITALNSIEGVYRKAVVLQAQNMSEA
HKLLQKLLVHCQKLKNTMEMVISVLLSVAELYWRSSSPTIALPMLLQALALSKEYRLQYLASETVLNLAFA
QLILGIPEQALSLLHMAIEPILADGAILDKGRAMFLVAKCQVASAASYDQPKKAEALEAAIENLNEAKNY
FAKVDCERIRDVVYFQARLYHTLGKTQERNRCAMLFRQLHQELPSHGVPLINHL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	84.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_057321
Locus ID:	51433
UniProt ID:	Q9UJX4
RefSeq Size:	2625
Cytogenetics:	12q24.31
RefSeq ORF:	2265
Synonyms:	APC5
Summary:	This gene encodes a tetratricopeptide repeat-containing component of the anaphase promoting complex/cyclosome (APC/C), a large E3 ubiquitin ligase that controls cell cycle progression by targeting a number of cell cycle regulators such as B-type cyclins for 26S proteasome-mediated degradation through ubiquitination. The encoded protein is required for the proper ubiquitination function of APC/C and for the interaction of APC/C with transcription coactivators. It also interacts with polyA binding protein and represses internal ribosome entry site-mediated translation. Multiple transcript variants encoding different isoforms have been found for this gene. These differences cause translation initiation at a downstream AUG and result in a shorter protein (isoform b), compared to isoform a. [provided by RefSeq, Nov 2008]
Protein Families:	Druggable Genome
Protein Pathways:	Cell cycle, Oocyte meiosis, Progesterone-mediated oocyte maturation, Ubiquitin mediated proteolysis

Product images:



Coomassie blue staining of purified ANAPC5 protein (Cat# [TP303030]). The protein was produced from HEK293T cells transfected with ANAPC5 cDNA clone (Cat# [RC203030]) using MegaTran 2.0 (Cat# [TT210002]).