

Product datasheet for **TP303079L**

PALMD (NM_017734) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human palmdelphin (PALMD), 1 mg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC203079 protein sequence
Red=Cloning site **Green**=Tags(s)

MEEAELVKGRQLQAITDKRKIQEEISQRLKIEEDKLKHQHLKKKALREKWLLDGISSGKEQEEMKKQNQQ
DQHQIQVLEQSILRLEKEIQDLEKAEQLISTKEEAILKKLSIERTTEDIIRSVKVEREERAEESEDY
ANIPDLPKSYIPSLRKEINEEKEDDEQNRKALYAMEIKVEKDLKTGSTVLSSIPLPSDDFKGTGKIVY
DDGQKSVYAVSSNHSAAANGTDGLAPVEVEELLRQASERNKSPTEYHEPVYANPFYRPTTPQRETVP
PNFQERIKIKTNGLGIGVNESIHNMGNLSEERGNNFNHISPIPPVPHPRSVIQQAEEKLHTPQKRLMTP
WEESNVMQDKDAPSPKPRLSPRETIFGKSEHQNSSPTCQEDEEDVRYNIVHSLPPDINDTEPVTMIFMGY
QQAEDSEEDKKFLTGYDGIHAELVVIDDEEEDEGEAEKPSYHPIAPHSQVYQPAKPTPLPRKRSEASP
HENTNHKSPHKNSISLKEQEESLGSPVHSPFDAQTTGDGTEDPSLTALRMRMAKLGKKVI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 62.6 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.


Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: NP_060204
 Locus ID: 54873
 UniProt ID: Q9NP74
 RefSeq Size: 2601
 Cytogenetics: 1p21.2
 RefSeq ORF: 1653
 Synonyms: C1orf11; PALML

Product images:


Coomassie blue staining of purified PALMD protein (Cat# [TP303079]). The protein was produced from HEK293T cells transfected with PALMD cDNA clone (Cat# [RC203079]) using MegaTran 2.0 (Cat# [TT210002]).