

## Product datasheet for **TP303086L**

### **GST3 (GSTP1) (NM\_000852) Human Recombinant Protein**

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human glutathione S-transferase pi 1 (GSTP1), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203086 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MPPYTVVYFPVRGRCAALRMLLADQGQSWKEEVTVETWQEGSLKASCLYGQLPKFQDGDLTLYQSNTIL RHLGRTLGLYGKDQQAALVDMVNDGVEDLRCKYVSLIYTNYEAGKDDYVKALPGQLKPFETLLSQNQGG KTFIVGDQISFADYNLLDLLLIHEVLAPGCLDAFPLLSAYVGRLSARPKLKAFLASPEYVNLPIGNGKQ  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	23.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_000843</a></u>
Locus ID:	2950
UniProt ID:	<u><a href="#">P09211</a></u> , <u><a href="#">V9HWE9</a></u>
RefSeq Size:	986



[View online »](#)

**Cytogenetics:** 11q13.2

**RefSeq ORF:** 630


**Synonyms:** DFN7; FAEES3; GST3; GSTP; HEL-S-22; PI

**Summary:** Glutathione S-transferases (GSTs) are a family of enzymes that play an important role in detoxification by catalyzing the conjugation of many hydrophobic and electrophilic compounds with reduced glutathione. Based on their biochemical, immunologic, and structural properties, the soluble GSTs are categorized into 4 main classes: alpha, mu, pi, and theta. This GST family member is a polymorphic gene encoding active, functionally different GSTP1 variant proteins that are thought to function in xenobiotic metabolism and play a role in susceptibility to cancer, and other diseases. [provided by RefSeq, Jul 2008]

**Protein Families:** Druggable Genome

**Protein Pathways:** Drug metabolism - cytochrome P450, Glutathione metabolism, Metabolism of xenobiotics by cytochrome P450, Pathways in cancer, Prostate cancer

### Product images:



Coomassie blue staining of purified GSTP1 protein (Cat# [TP303086]). The protein was produced from HEK293T cells transfected with GSTP1 cDNA clone (Cat# [RC203086]) using MegaTran 2.0 (Cat# [TT210002]).