

Product datasheet for TP303114L

ELK3 (NM_005230) Human Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant protein of human ELK3, ETS-domain protein (SRF accessory protein 2) (ELK3), 1 mg |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | >RC203114 protein sequence Red=Cloning site Green=Tags(s) |

MESAITLWQFLLQLLLDQKHEHLICWTSNDGEFKLLKAEVAKLWGLRKNKTNMNYDKLSRALRYYYDKN
IHKVIGQKFVYKFSVFEILKMDPHAVEISRESLLLQSDCKASPEGREAHKHGLAALRSTSRNEYIHS
GLYSSFTINSLQNPPDAFKAIKTEKLEPPEDSPPVEEVRTVIRFVTNKTDKHVTRPVVSLPSTSEAAAA
SAFLASSVSAKISSLMLPNAASISSASPFSSRSPSLSPNSPLPSEHRS LFLEAACHSDSDSLEPLNLSSGS
KTKSPSLPPKAKKPKGLEISAPPLVLSGTDIGSIALNSPALPSGSLTPAFFTAQTPNGLLLTPSPLSSI
HFWSSLPVAPLSPARLQGPSTLFQFPTLLNGHMPVPIPSLDRAASPVLLSSNSQKS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

| | |
|----------------|--|
| Tag: | C-Myc/DDK |
| Predicted MW: | 44.1 kDa |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol |
| Preparation: | Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps. |
| Note: | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage: | Store at -80°C. |
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| RefSeq: | <u>NP_005221</u> |




[View online »](#)

Locus ID: 2004
UniProt ID: [P41970](#), [A0A024RBE2](#)
RefSeq Size: 2180
Cytogenetics: 12q23.1
RefSeq ORF: 1221
Synonyms: ERP; NET; SAP-2; SAP2

Summary: This gene encodes a member of the ETS-domain transcription factor family and the ternary complex factor (TCF) subfamily. Proteins in this subfamily regulate transcription when recruited by serum response factor to bind to serum response elements. This protein is activated by signal-induced phosphorylation; studies in rodents suggest that it is a transcriptional inhibitor in the absence of Ras, but activates transcription when Ras is present. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Jan 2015]

Protein Families: Transcription Factors

Product images:



Coomassie blue staining of purified ELK3 protein (Cat# [TP303114]). The protein was produced from HEK293T cells transfected with ELK3 cDNA clone (Cat# [RC203114]) using MegaTran 2.0 (Cat# [TT210002]).