

Product datasheet for TP303128

Tetranectin (CLEC3B) (NM_003278) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human C-type lectin domain family 3, member B (CLEC3B), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203128 protein sequence Red =Cloning site Green =Tags(s)
	MELWGAYLLLCLFSLTQVTEPPTQKPKKIVNAKDVVNTKMFEEELKSRLDTLAQEVALLEKQQALQTV CLKGTVHMKCFLAFTQTKTFHEASEDCISRGGTLSTPQTGSENDALYEYLRQSVGNEAEIWLGLNDMAA EGTWVDMTGARIAYKNWETEITAQPDGGKTENCAVLGAANGKWFDRKCRDQLPYICQFGIV
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	22.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_003269
Locus ID:	7123
UniProt ID:	P05452
RefSeq Size:	886



[View online »](#)

Cytogenetics: 3p21.31


RefSeq ORF: 606

Synonyms: TN; TNA

Summary: Tetranectin binds to plasminogen and to isolated kringle 4. May be involved in the packaging of molecules destined for exocytosis.[UniProtKB/Swiss-Prot Function]

Protein Families: Secreted Protein

Product images:



Coomassie blue staining of purified CLEC3B protein (Cat# TP303128). The protein was produced from HEK293T cells transfected with CLEC3B cDNA clone (Cat# [RC203128]) using MegaTran 2.0 (Cat# [TT210002]).