

Product datasheet for TP303174L

Glutathione Synthetase (GSS) (NM_000178) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human glutathione synthetase (GSS), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203174 protein sequence Red=Cloning site Green=Tags(s)

MATNWGSELLQDKQQLLEELARQAVDRALAEVLLRSTSSEVVSYPFTLFPVSLVPSALLEQAYAVQM
DFNLLVDAVSQNAAFLEQTLSSSTIKQDDFTARLFDIHKQVLKEGIAQTVFLGLNRSYDMFQRSADGSPAL
KQIEINTISASFGGLASRTPAVHRHVLSVLSKTKEAGKILSNPSKGLALGIAKAWELYGSPNALVLLIA
QEKERNIFDQRAIENELLARNIHVIRRTFEDISEKGLDQDRRLFVDGQEIADVYFRDGYMQRQYSLQNW
EARLLERSHAAKCPDIATQLAGTKKVQQLSRPGMLEMLLPQPEAVARLRATFAGLYSLDVGEEGDQA
IAEALAAPSRFVLKPQREGGNNLYGEEMVQALKQLKDSEERASYILMEKIEPEPFENCLLRPGSPARVV
QCISELGFVYVRQEKTLVMNKHVGHLLRRTKAIEHADGGVAAGVAVLDPYPV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	52.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_000169



[View online »](#)

Locus ID:	2937
UniProt ID:	P48637 , V9HWJ1
RefSeq Size:	1918
Cytogenetics:	20q11.22
RefSeq ORF:	1422
Synonyms:	GSHS; HEL-S-64p; HEL-S-88n
Summary:	Glutathione is important for a variety of biological functions, including protection of cells from oxidative damage by free radicals, detoxification of xenobiotics, and membrane transport. The protein encoded by this gene functions as a homodimer to catalyze the second step of glutathione biosynthesis, which is the ATP-dependent conversion of gamma-L-glutamyl-L-cysteine to glutathione. Defects in this gene are a cause of glutathione synthetase deficiency. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome
Protein Pathways:	Glutathione metabolism, Metabolic pathways

Product images:



Coomassie blue staining of purified GSS protein (Cat# [TP303174]). The protein was produced from HEK293T cells transfected with GSS cDNA clone (Cat# [RC203174]) using MegaTran 2.0 (Cat# [TT210002]).