

Product datasheet for **TP303202L**

EEF1D (NM_001960) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human eukaryotic translation elongation factor 1 delta (guanine nucleotide exchange protein) (EEF1D), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203202 protein sequence Red =Cloning site Green =Tags(s) MRSKGASCTLETVWEDKHKYEEAERRFYEHEATQAAASAQQLPAEGPAMNGPGQDDPEDADEAEAPDGGS RRDPRKSQDSRKPLQKKRKRSPKSGLGPADLALLGLSAERVWLDKSLFDQAESSYRQKLADVAAQAAWPP ALAPWGLCTHGNQVACHHVTWGIWVNKSSFDQAERAFVEWSQALLAPEGSRROGTPNTGQQVAVPDLAH QPSPVNGQPPLGSLQALVREVWLEKPRYDAAERGFYEALFDGHPGKVRQLERAGLAEGARRRRDRRG RNILGNKRAGLRRADGEAPSALPYCYFLQKDAEAPWLSKPAYDSAECRHAAEALRVAWCLEAASLSHRP GPRSGLSVSSLRPNRKMATNFLAHEKIWFDFKFKYDDAERRFYEQMNGPVAGASRQENGASVILRDIARAR ENIQKSLAGSSGPGASSGTSGDHGELVRIASLEVENQSLRGVVQELQQAIKLEARLNVLEKSSPGHRA TAPQTQHVSMPMRQVEPPAKKPATPAEDDEDDEDIDLFGSDNEEDKEAAQLREERLRQYAEKKAKKPALVA KSSILLDVKPWDEDETDMAQLEACVRSIQLDGLVWGASKLVPVGYGIRKLQIQCVVEDDKVGTDLLEEEIT KFEHVQSVDAIAFNKI TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	30.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Bioactivity:	Biolayer interferometry (BLI) assay (PMID: 26624286)
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.




[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001951</u>
Locus ID:	1936
UniProt ID:	<u>P29692</u>
RefSeq Size:	1276
Cytogenetics:	8q24.3
RefSeq ORF:	1944
Synonyms:	EF-1D; EF1D; FP1047
Summary:	This gene encodes a subunit of the elongation factor-1 complex, which is responsible for the enzymatic delivery of aminoacyl tRNAs to the ribosome. This subunit, delta, functions as guanine nucleotide exchange factor. It is reported that following HIV-1 infection, this subunit interacts with HIV-1 Tat. This interaction results in repression of translation of host cell proteins and enhanced translation of viral proteins. Several alternatively spliced transcript variants encoding multiple isoforms have been found for this gene. Related pseudogenes have been defined on chromosomes 1, 6, 7, 9, 11, 13, 17, 19.[provided by RefSeq, Aug 2010]

Product images:



Biolayer interferometry assays (BLI) of a biotinylated HIV RTp66/p51-loaded or integrase-loaded biosensor with 60 nM of purified eEF1A (OriGene [TP309697]), eEF1B (OriGene [TP321264]), eEF1D (OriGene [TP303202]), eEF1G (OriGene [TP319578]) or BSA as analytes. Figure cited from PLoS Pathog, PMID: 26624286



Coomassie blue staining of purified EEF1D protein (Cat# [TP303202]). The protein was produced from HEK293T cells transfected with EEF1D cDNA clone (Cat# [RC203202]) using MegaTran 2.0 (Cat# [TT210002]).