

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TP303240M

MHF2 (CENPX) (NM_144998) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins	
Description:	Recombinant protein of human stimulated by retinoic acid 13 homolog (mouse) (STRA13), 100 µg	
Species:	Human	
Expression Host:	HEK293T	
Expression cDNA Clone or AA Sequence:	>RC203240 protein sequence Red=Cloning site Green=Tags(s)	
	MEGAGAGSGFRKELVSRLLHLHFKDDKTKEAAVRGVRQAQAEDALRADVDQLEKVLPQLLLDF	
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV	
Tag:	C-Myc/DDK	
Predicted MW:	6.8 kDa	
Concentration:	>0.05 µg/µL as determined by microplate BCA method	
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining	
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol	
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.	
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.	
Storage:	Store at -80°C.	
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.	
RefSeq:	<u>NP 659435</u>	
Locus ID:	201254	
UniProt ID:	<u>A8MT69</u>	
RefSeq Size:	875	
Cytogenetics:	17q25.3	



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2025 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	MHF2 (CENPX) (NM_144998) Human Recombinant Protein – TP303240M
RefSeq ORF:	189
Synonyms:	CENP-X; D9; FAAP10; MHF2; STRA13
Summary:	DNA-binding component of the Fanconi anemia (FA) core complex. Required for the normal activation of the FA pathway, leading to monoubiquitination of the FANCI-FANCD2 complex in response to DNA damage, cellular resistance to DNA cross-linking drugs, and prevention of chromosomal breakage (PubMed:20347428, PubMed:20347429). In complex with CENPS (MHF heterodimer), crucial cofactor for FANCM in both binding and ATP-dependent remodeling of DNA. Stabilizes FANCM. In complex with CENPS and FANCM (but not other FANC proteins), rapidly recruited to blocked forks and promotes gene conversion at blocked replication forks (PubMed:20347428, PubMed:20347429). In complex with CENPS, CENPT and CENPW (CENP-T-W-S-X heterotetramer), involved in the formation of a functional kinetochore outer plate, which is essential for kinetochore-microtubule attachment and faithful mitotic progression (PubMed:19620631). As a component of MHF and CENP-T-W-S-X complexes, binds DNA and bends it to form a nucleosome-like structure (PubMed:20347428, PubMed:20347429). DNA-binding function is fulfilled in the presence of CENPS, with the following preference for DNA substates: Holliday junction > double-stranded > splay arm > single-stranded. Does not bind DNA on its own (PubMed:20347429).[UniProtKB/Swiss-Prot Function]

Product images:

116	_	
66	_	
45	_	
35	_	
25	_	
18 14	=	-

Coomassie blue staining of purified CENPX protein (Cat# [TP303240]). The protein was produced from HEK293T cells transfected with CENPX cDNA clone (Cat# [RC203240]) using MegaTran 2.0 (Cat# [TT210002]).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2025 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US