

## Product datasheet for **TP303250M**

### Phospholipase A2 IIA (PLA2G2A) (NM\_000300) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human phospholipase A2, group IIA (platelets, synovial fluid) (PLA2G2A), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203250 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MKTLVVAVIMIFGLLQAHGNLVNFHRMIKLTGKEAALSYGFYGCCHGCVGGGRGSPKDATDRCCVTHDCC YKRLEKRGCGTKFLSYKFSNSGRITCAKQDSCRSQLCECDKAAATCFARNKTTYNKKYQYYSNKHCRGS TPRC
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	13.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_000291</a>
Locus ID:	5320
UniProt ID:	<a href="#">P14555</a> , <a href="#">A0A024RA96</a>



[View online »](#)

RefSeq Size: 1017

Cytogenetics: 1p36.13

RefSeq ORF: 432


Synonyms: MOM1; PLA2; PLA2B; PLA2L; PLA2S; PLAS1; sPLA2

**Summary:** The protein encoded by this gene is a member of the phospholipase A2 family (PLA2). PLA2s constitute a diverse family of enzymes with respect to sequence, function, localization, and divalent cation requirements. This gene product belongs to group II, which contains secreted form of PLA2, an extracellular enzyme that has a low molecular mass and requires calcium ions for catalysis. It catalyzes the hydrolysis of the sn-2 fatty acid acyl ester bond of phosphoglycerides, releasing free fatty acids and lysophospholipids, and thought to participate in the regulation of the phospholipid metabolism in biomembranes. Several alternatively spliced transcript variants with different 5' UTRs have been found for this gene. [provided by RefSeq, Sep 2009]

**Protein Families:** Druggable Genome, Transmembrane

**Protein Pathways:** alpha-Linolenic acid metabolism, Arachidonic acid metabolism, Ether lipid metabolism, Fc epsilon RI signaling pathway, Glycerophospholipid metabolism, GnRH signaling pathway, Linoleic acid metabolism, Long-term depression, MAPK signaling pathway, Metabolic pathways, Vascular smooth muscle contraction, VEGF signaling pathway

### Product images:



Coomassie blue staining of purified PLA2G2A protein (Cat# [TP303250]). The protein was produced from HEK293T cells transfected with PLA2G2A cDNA clone (Cat# [RC203250]) using MegaTran 2.0 (Cat# [TT210002]).