

Product datasheet for **TP303257L**

MT1H (NM_005951) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human metallothionein 1H (MT1H), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203257 protein sequence Red =Cloning site Green =Tags(s)
	MDPNCSC EAGGSCACAGSCKCKKCKCTSCCKSCCSCCPLGCAKCAQGCICKGASEKCSCCA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	5.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_005942
Locus ID:	4496
UniProt ID:	P80294 , A0A140VJP7
RefSeq Size:	397
Cytogenetics:	16q13



[View online »](#)

RefSeq ORF:	183
Synonyms:	MT-0; MT-1H; MT-IH; MT1
Summary:	Metallothioneins have a high content of cysteine residues that bind various heavy metals; these proteins are transcriptionally regulated by both heavy metals and glucocorticoids. [UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified MT1H protein (Cat# [TP303257]). The protein was produced from HEK293T cells transfected with MT1H cDNA clone (Cat# [RC203257]) using MegaTran 2.0 (Cat# [TT210002]).