

Product datasheet for **TP303264M**

NIPSNAP3A (NM_015469) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human nipsnap homolog 3A (C. elegans) (NIPSNAP3A), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203264 protein sequence Red =Cloning site Green =Tags(s)
	<p>MLVLRSA LTRALASRTLAPQMCSSFATGPRQYDGIFYEFRSYLKP SKMNEFLENFEKNAHLRTAHSELV GYWSVEFGGRMNTVFHIWKYDNFAHRTEVRKALAKDKEWQEQLIPNLALIDKQSEITYLVPWCKLEKP PKEGVYELATFQMKPGGPALWGDAFKRAVHAHVNLGYTKLVGVFHTTEYGALNRVHVLWWNESADSRRAA GR HKSHEDPRVVAAVRESVNYLVSQQNMLLIPTSFSP LK</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	28.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_056284</u>
Locus ID:	25934



[View online »](#)

UniProt ID:	<u>Q9UFN0</u>
RefSeq Size:	1664
Cytogenetics:	9q31.1
RefSeq ORF:	741
Synonyms:	HSPC299; NIPSNAP4; TASSC
Summary:	NIPSNAP3A belongs to a family of proteins with putative roles in vesicular transport (Buechler et al., 2004 [PubMed 15177564]).[supplied by OMIM, Mar 2008]

Product images:



Coomassie blue staining of purified NIPSNAP3A protein (Cat# [TP303264]). The protein was produced from HEK293T cells transfected with NIPSNAP3A cDNA clone (Cat# [RC203264]) using MegaTran 2.0 (Cat# [TT210002]).