

Product datasheet for TP303286

C6orf48 (NM_001040438) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 6 open reading frame 48 (C6orf48), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203286 protein sequence Red=Cloning site Green=Tags(s) MERSFVWLSCLDSDSCNLTFRLGVEVESHACSPSLLWNLLTQYLPPGAGHILRTYNFPVLSCVSSCHLIGG KMPEN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	8.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001035528
Locus ID:	50854
RefSeq Size:	751
Cytogenetics:	6p21.33




[View online »](#)

RefSeq ORF: 225

Synonyms: D6S57; G8

Product images:



Coomassie blue staining of purified C6orf48 protein (Cat# TP303286). The protein was produced from HEK293T cells transfected with C6orf48 cDNA clone (Cat# [RC203286]) using MegaTran 2.0 (Cat# [TT210002]).