

Product datasheet for TP303288

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

MAGEA3 (NM_005362) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human melanoma antigen family A, 3 (MAGEA3), 20 μg

Species: Human
Expression Host: HEK293T

Expression cDNA Clone >RC203288 representing NM_005362 **or AA Sequence:** Red=Cloning site Green=Tags(s)

MPLEQRSQHCKPEEGLEARGEALGLVGAQAPATEEQEAASSSSTLVEVTLGEVPAAESPDPPQSPQGASS LPTTMNYPLWSQSYEDSSNQEEEGPSTFPDLESEFQAALSRKVAELVHFLLLKYRAREPVTKAEMLGSVV GNWQYFFPVIFSKASSSLQLVFGIELMEVDPIGHLYIFATCLGLSYDGLLGDNQIMPKAGLLIIVLAIIA REGDCAPEEKIWEELSVLEVFEGREDSILGDPKKLLTQHFVQENYLEYRQVPGSDPACYEFLWGPRALVE

TSYVKVLHHMVKISGGPHISYPPLHEWVLREGEE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 34.6 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by

conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 005353

Locus ID: 4102






UniProt ID: P43357 1753 RefSeq Size: Cytogenetics: Xq28 942 RefSeq ORF:

Synonyms: CT1.3; HIP8; HYPD; MAGE3; MAGEA6

Summary: This gene is a member of the MAGEA gene family. The members of this family encode

proteins with 50 to 80% sequence identity to each other. The promoters and first exons of the MAGEA genes show considerable variability, suggesting that the existence of this gene family enables the same function to be expressed under different transcriptional controls. The MAGEA genes are clustered at chromosomal location Xq28. They have been implicated in some hereditary disorders, such as dyskeratosis congenita. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified MAGEA3 protein (Cat# TP303288). The protein was produced from HEK293T cells transfected with MAGEA3 cDNA clone (Cat# [RC203288]) using MegaTran 2.0 (Cat# [TT210002]).