

Product datasheet for **TP303293**

HSD17B2 (NM_002153) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hydroxysteroid (17-beta) dehydrogenase 2 (HSD17B2), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203293 protein sequence Red =Cloning site Green =Tags(s)

MSTFFSDTAWICLAVPTVLCGTVFCKYKKSSGQLWSWMVCLAGLCAVCLLILSPFWGLILFSVSCFLMYT
YLSGQELLPVDQKAVLVTGGDCGLGHALCKYLDELGFTVFAGVLNENGPGEELRRTCSPRLSVLQMDIT
KPVQIKDAYSKVAAMLQDRGLWAVINNAGVLGFPTDGELLMTDYKQCMAVNFFGTVEVTKTFLPLLRKS
KGRLVNVSSMGGGAPMERLASYGSSKAAVTMFSSVMRLELSKWGIKVASIQPGGFLTNIAGTSDKWEKLE
KDILDHLP AEVQEDYGQDYILAQRNLLLLINSLASKDFSPVLRDIQHAILAKSPFAYYTPGKGAYLWICL
AHYLPIGIYDYFAKRHFGQDKPMPRALRMPNYKKKAT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	42.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_002144
Locus ID:	3294




[View online »](#)

UniProt ID:	P37059
RefSeq Size:	1451
Cytogenetics:	16q23.3
RefSeq ORF:	1161
Synonyms:	EDH17B2; HSD17; SDR9C2
Summary:	Capable of catalyzing the interconversion of testosterone and androstenedione, as well as estradiol and estrone. Also has 20-alpha-HSD activity. Uses NADH while EDH17B3 uses NADPH.[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Androgen and estrogen metabolism, Metabolic pathways

Product images:



HSD17B2 enzymatic activity with 75µM β-estradiol as a substrate, measured by NADH production (indicated by luminescence).



Coomassie blue staining of purified HSD17B2 protein (Cat# TP303293). The protein was produced from HEK293T cells transfected with HSD17B2 cDNA clone (Cat# [RC203293]) using MegaTran 2.0 (Cat# [TT210002]).