

## **Product datasheet for TP303301L**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## C1D (NM\_173177) Human Recombinant Protein

**Product data:** 

**Product Type:** Recombinant Proteins

**Description:** Recombinant protein of human C1D nuclear receptor co-repressor (C1D), transcript variant 2,

1 mg

Species: Human
Expression Host: HEK293T

**Expression cDNA Clone** >RC203301 protein sequence or AA Sequence: Red=Cloning site Green=Tags(s)

MAGEEINEDYPVEIHEYLSAFENSIGAVDEMLKTMMSVSRNELLQKLDPLEQAKVDLVSAYTLNSMFWVY LATQGVNPKEHPVKQELERIRVYMNRVKEITDKKKAGKLDRGAASRFVKNALWEPKSKNASKVANKGKSK

S

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV** 

Tag: C-Myc/DDK
Predicted MW: 15.8 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining

**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

**Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by

conventional chromatography steps.

**Note:** For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

**RefSeq:** NP 775269

Locus ID: 10438
UniProt ID: <u>Q13901</u>





RefSeq Size: 1194

Cytogenetics: 2p14 RefSeq ORF: 423

Synonyms: hC1D; LRP1; Rrp47; SUN-CoR; SUNCOR


**Summary:** The protein encoded by this gene is a DNA binding and apoptosis-inducing protein and is

localized in the nucleus. It is also a Rac3-interacting protein which acts as a corepressor for the thyroid hormone receptor. This protein is thought to regulate TRAX/Translin complex formation. Alternate splicing results in multiple transcript variants that encode the same protein. Multiple pseudogenes of this gene are found on chromosome 10.[provided by

RefSeq, Jun 2010]

Protein Families: Druggable Genome
Protein Pathways: RNA degradation

## **Product images:**



Coomassie blue staining of purified C1D protein (Cat# [TP303301]). The protein was produced from HEK293T cells transfected with C1D cDNA clone (Cat# [RC203301]) using MegaTran 2.0 (Cat# [TT210002]).