

Product datasheet for **TP303390**

ZNF447 (ZSCAN18) (NM_023926) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Purified recombinant protein of Homo sapiens zinc finger and SCAN domain containing 18 (ZSCAN18), 20 µg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC203390 protein sequence
Red=Cloning site **Green**=Tags(s)

MLPLEKAFASPRSSPAPPDLTPGSAAGVQEEPETIPERTPADLEFSRLRFREFVYQEAAGPHQTLARL
HELCRQWLMPEARSKEQMLELLVLEQFLGILPDKVRPWVVAQYPESCCKAASLVEGLADVLEEPGMMLLGS
PAGSSSILSDGVYERHMDPLLLPGELASPSQALGAGEIPAPSETPWLSPDPLFLEQRRVREAKTEEDGPA
NTEQKLKSPEDPQHLGEWGHLDPAEENLKSyrKLLLWGyQLSQPDAASRLDTEELRLVERDPQGSSSLPE
GGRRQESAGCACEEAAPAGVLPPELPTTEAPPGDALADPPSGTTEEEEEEQPGKAPDPQDPQDAESDSATGSQ
RQSVIQQPAPDRGTAKLGTKRPHPEDGDEQSLEGVSSSGDSAGLEAGQGPGADEPGLSRGKPYACGECGE
AFAWLSHLMESHSSHGGRKRYACQGCWKTFFHSLALAEHQKTHEKEKSYALGGARGPQPSTREAQAGARA
GGPPESVEGEAPPAPPEAQR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 54.6 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.


Storage: Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_076415
Locus ID:	65982
UniProt ID:	Q8TBC5 , A0A024R4T0
RefSeq Size:	2960
Cytogenetics:	19q13.43
RefSeq ORF:	1530
Synonyms:	ZNF447
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Protein Families:	Transcription Factors

Product images:



Coomassie blue staining of purified ZSCAN18 protein (Cat# TP303390). The protein was produced from HEK293T cells transfected with ZSCAN18 cDNA clone (Cat# [RC203390]) using MegaTran 2.0 (Cat# [TT210002]).