

Product datasheet for TP303403

UBTD1 (NM_024954) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ubiquitin domain containing 1 (UBTD1), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203403 representing NM_024954 Red=Cloning site Green=Tags(s) MGNCVGRQRRERPAAPGHPKRKRAGRNEPLKKERLKWKSDYPMTDGGQLRSKRDEFWDTAPAFEGRKEIWDA LKAAAYAAEANDHELAQAILDGASITLPHGTLCECYDELGNRYQLPIYCLSPVNLLEHTEESLEPPE PPPSVRREFPLKVRLSTGKDVRLSASLPDVTGQLKRQLHAQEGIEPSWQRWFFSGKLLTDRTRLQETKIQ KDFVIQVIINQPPPPQD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	25.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_079230
Locus ID:	80019
UniProt ID:	Q9HAC8



[View online »](#)

RefSeq Size: 1575

Cytogenetics: 10q24.1-q24.2

RefSeq ORF: 681

Summary: The degradation of many proteins is carried out by the ubiquitin pathway in which proteins are targeted for degradation by covalent conjugation of the polypeptide ubiquitin. This gene encodes a protein that belongs to the ubiquitin family of proteins. The encoded protein is thought to regulate E2 ubiquitin conjugating enzymes belonging to the UBE2D family. [provided by RefSeq, Mar 2014]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified UBTD1 protein (Cat# TP303403). The protein was produced from HEK293T cells transfected with UBTD1 cDNA clone (Cat# [RC203403]) using MegaTran 2.0 (Cat# [TT210002]).