

## Product datasheet for TP303412L

### RPL13 (NM\_033251) Human Recombinant Protein

#### Product data:

|                                       |   |
|---------------------------------------|---|
| Product Type:                         | Recombinant Proteins  |
| Description:                          | Recombinant protein of human ribosomal protein L13 (RPL13), transcript variant 2, 1 mg  |
| Species:                              | Human   |
| Expression Host:                      | HEK293T   |
| Expression cDNA Clone or AA Sequence: | >RC203412 protein sequence<br><b>Red</b> =Cloning site <b>Green</b> =Tags(s)  |
|                                       | <br>MAPSRNGMVLKPHFHKDWQRRVATWFNQPARKIRRRKARQAKARRIAPRPASGPIRPIVRCPTVRYHTKV<br>RAGRGFSLEELRVAGIHKKVARTIGISVDPRRRNKSTESLQANVQLKEYRSKLILFPRKPSAPKKGDSS<br>AEELKLATQLTGPVMPVVRNVYKKEKARVITEEEKNFKAFASLRMARANARLFGIRAKRAKEAAEQDVEKK<br>K<br><br><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b> |
| Tag:                                  | C-Myc/DDK   |
| Predicted MW:                         | 24.1 kDa  |
| Concentration:                        | >0.05 µg/µL as determined by microplate BCA method  |
| Purity:                               | > 80% as determined by SDS-PAGE and Coomassie blue staining   |
| Buffer:                               | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol  |
| Preparation:                          | Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.  |
| Note:                                 | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.  |
| Storage:                              | Store at -80°C.   |
| Stability:                            | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.   |
| RefSeq:                               | <a href="#">NP_150254</a>   |
| Locus ID:                             | 6137  |
| UniProt ID:                           | <a href="#">P26373</a> , <a href="#">A8K4C8</a>   |



[View online »](#)

RefSeq Size: 4672

Cytogenetics: 17p11.2

RefSeq ORF: 633


Synonyms: BBC1; D16S44E; D16S444E; L13; SEMDIST

**Summary:** Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L13E family of ribosomal proteins. It is located in the cytoplasm. This gene is expressed at significantly higher levels in benign breast lesions than in breast carcinomas. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. [provided by RefSeq, Jul 2011]

**Protein Families:** Druggable Genome

**Protein Pathways:** Ribosome

### Product images:



Coomassie blue staining of purified RPL13 protein (Cat# [TP303412]). The protein was produced from HEK293T cells transfected with RPL13 cDNA clone (Cat# [RC203412]) using MegaTran 2.0 (Cat# [TT210002]).