

Product datasheet for **TP303441L**

COPS4 (NM_016129) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human COP9 constitutive photomorphogenic homolog subunit 4 (Arabidopsis) (COPS4), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203441 protein sequence Red =Cloning site Green =Tags(s)

MAAAVRQDLAQLMNSSGSHKDLGKYRQILEKAIQLSGAEQLEALKAFVEAMVNVSLVISRQLLTD
THLPNLPDSTAKEIYHFTLEKIQRVISFEEQVASIRQHLSIYEKEEDWRNAAQVLVGIPL
VDYKLETYLKIARLYLEDDDPVQAEAYINRASLLQNESTNEQLQIHYKVCYARVLDYRRK
LSYKTIVHESERLEALKHALHCTILASAGQQRSMRLATLFKDERCQQLAAYGILEKMYLDRI
IRGNQLQEFAAMLMPHQKATTADGSSILDRAVIEHNLLSASKLYNNITFEELGALLEIPA
AKAEKIASQMITEGRMNGFIDQIDGIVHFETREALPTWQKIQSLCFQVNNLLEKISQTAP
EWTAQAMEAQMAQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	46.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_057213</u>




[View online »](#)

Locus ID: 51138
UniProt ID: [Q9BT78](#)
RefSeq Size: 1765
Cytogenetics: 4q21.22
RefSeq ORF: 1218
Synonyms: CSN4; SGN4

Summary: This gene encodes one of eight subunits composing COP9 signalosome, a highly conserved protein complex that functions as an important regulator in multiple signaling pathways. The structure and function of COP9 signalosome is similar to that of the 19S regulatory particle of 26S proteasome. COP9 signalosome has been shown to interact with SCF-type E3 ubiquitin ligases and act as a positive regulator of E3 ubiquitin ligases. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2012]

Product images:



Coomassie blue staining of purified COPS4 protein (Cat# [TP303441]). The protein was produced from HEK293T cells transfected with COPS4 cDNA clone (Cat# [RC203441]) using MegaTran 2.0 (Cat# [TT210002]).