

## Product datasheet for TP303490

### ATIC (NM\_004044) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human 5-aminoimidazole-4-carboxamide ribonucleotide formyltransferase/IMP cyclohydrolase (ATIC), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203490 protein sequence Red=Cloning site Green=Tags(s)

MAPGQLALFSVSDKTGLVEFARNLTALGLNLVASGGTAKALRDAGLAVRDVSELTFPPEMLGGRVKTLLHP  
AVHAGILARNIPEDNADMARLDFNLIRVVACNLYPFKTVASPGVTVEEAVEQIDIGGVTLLRAAAKNHA  
RVTVVCEPEDYVWVSTEMQSSKSDTSLETRRQLALKAFTHTAQYDEAISDYFRKQYKSGVSMPLRYGM  
NPHQTPAQLYTLQPKLPITVLNGAPGFINLDCDALNAWQLVKELKEALGIPAAASFKHVSPAGAAVGIPLS  
EDEAKVCMVYDLYKTLTPISAAYARARGADRMSSFGDFVALSDVCDVPTAKIISREVS DGIIAPGYEEEE  
LTILSKKKNGNYCVLQMDQSYKPDENEVRTLFGHLHSQKRNNGWVDKSLFSNVVTKNKDLPESALRDLIV  
ATIAVKYTQNSVICYAKNGQVIGIGAGQQSRIHCTRLAGDKANYWWLRHHPQVLSMKFKTGVKRAEISNA  
IDQYVTGTIGEDDELKWKALFEEVPELLTEAEKKEWVEKLTVEVSISSDAFFPFRDNVDRAKRSGVAYIA  
APSGSAADKVIEACDELGIIAHTNLRLFFFH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	64.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_004035</a>
<b>Locus ID:</b>	471
<b>UniProt ID:</b>	<a href="#">P31939</a> , <a href="#">V9HWH7</a>
<b>RefSeq Size:</b>	2094
<b>Cytogenetics:</b>	2q35
<b>RefSeq ORF:</b>	1776
<b>Synonyms:</b>	AICAR; AICARFT; HEL-S-70p; IMPCHASE; PURH
<b>Summary:</b>	This gene encodes a bifunctional protein that catalyzes the last two steps of the de novo purine biosynthetic pathway. The N-terminal domain has phosphoribosylaminoimidazolecarboxamide formyltransferase activity, and the C-terminal domain has IMP cyclohydrolase activity. A mutation in this gene results in AICA-ribosiduria. [provided by RefSeq, Sep 2009]
<b>Protein Families:</b>	Stem cell - Pluripotency
<b>Protein Pathways:</b>	Metabolic pathways, One carbon pool by folate, Purine metabolism

### Product images:



Coomassie blue staining of purified ATIC protein (Cat# TP303490). The protein was produced from HEK293T cells transfected with ATIC cDNA clone (Cat# [RC203490]) using MegaTran 2.0 (Cat# [TT210002]).