

Product datasheet for TP303508

Neuropilin 1 (NRP1) (NM_001024628) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human neuropilin 1 (NRP1), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203508 protein sequence Red =Cloning site Green =Tags(s)

MERGLPLLCAVLALVLPAGAFRNDKCGDTIKIESPGYLTSPGYPHSYHPSEKCEWLIQAPDPYQRIMIN
FNPHFDLEDRDCKYDYVEVFDGENENGHFRGKFCGKIAPPPVSSGPFLFIKFVSDYETHGAGFSIRYEI
FKRGPECSQNYTTPSGVIKSPGFPEKYPNSLECTYIVFAPKMSEIILEFESFDLEPDSNPPGGMFRCYDR
LEIWDGFPDVGPHIGRYCGQKTPGRIRSSSGILSMVFYDLSAIAKEGFSANYSVLQSSVSEDFKCMEALG
MESGEIHSQITASSQYSTNWSAERSRLNYPENGWTPGEDSYREWIQVDLGLLRFRVAVGTQGAISKETK
KKYYVKTYKIDVSSNGEDWITIKEGNKPVLFGNTNPTDWWAVFPKPLITRFVRIKPATWETGISMRFE
VYGCKITDYPCSGMLGMVSGGLISDSQITSSNQGDRNWMPENIRLVTSRSGWALPPAPHSYINEWLQIDLG
EEKIVRGIHQGGKHRENKVFMRKFKIGYSNNGSDWKMIMDDSKRKAKEFEGNNNYDTPELRTFPALSTR
FIRIYPERATHGGLGLRMELLGCEVEAPTAGPTTPNGNLVDECDDQANCHSGTGDDDFQLTGGTTVLATE
KPTVIDSTIQSGIK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	71.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_001019799](#)

Locus ID: 8829

UniProt ID: [O14786](#), [Q68DN3](#), [O14786-2](#), [C9JJV3](#)

RefSeq Size: 2478

Cytogenetics: 10p11.22

RefSeq ORF: 1932


Synonyms: BDCA4; CD304; NP1; NRP; VEGF165R

Summary: This gene encodes one of two neuropilins, which contain specific protein domains which allow them to participate in several different types of signaling pathways that control cell migration. Neuropilins contain a large N-terminal extracellular domain, made up of complement-binding, coagulation factor V/VIII, and meprin domains. These proteins also contains a short membrane-spanning domain and a small cytoplasmic domain. Neuropilins bind many ligands and various types of co-receptors; they affect cell survival, migration, and attraction. Some of the ligands and co-receptors bound by neuropilins are vascular endothelial growth factor (VEGF) and semaphorin family members. This protein has also been determined to act as a co-receptor for SARS-CoV-2 (which causes COVID-19) to infect host cells. [provided by RefSeq, Nov 2020]

Protein Families: Druggable Genome, Secreted Protein, Transmembrane

Protein Pathways: Axon guidance

Product images:



Coomassie blue staining of purified NRP1 protein (Cat# TP303508). The protein was produced from HEK293T cells transfected with NRP1 cDNA clone (Cat# [RC203508]) using MegaTran 2.0 (Cat# [TT210002]).