

Product datasheet for **TP303513L**

TMEM41A (NM_080652) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human transmembrane protein 41A (TMEM41A), 1 mg

Species: Human

Expression Host: HEK293T

**Expression cDNA Clone
or AA Sequence:** >RC203513 protein sequence
Red=Cloning site **Green**=Tags(s)

MRPLLGLLLVFAGCTFALYLLSTRLPRGRRLGSTEEAGGRSLWFPSDLAELRESEVLREYRKEHQAYVF
LLFCGAYLYKQGFAIPGSSFLNVLGALFGPWLGLLLCCVLTSVGATCCYLLSSIFGKQLVVSYPDKVA
LLQRKVEENRNSLFFLLFLRLFPMTPNWFLNLSAPILNIPIVQFFFSVLIGLIPYNFICVQTGSILSTL
TSLDALFSWDTVFKLLAIAMVALIPGTLIKKFSQKHLQLNETSTANHIHSRKDT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 29.5 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_542383](#)


Locus ID: 90407

UniProt ID: [Q96HV5](#)



[View online »](#)

RefSeq Size:	2840
Cytogenetics:	3q27.2
RefSeq ORF:	792
Synonyms:	2900010K02Rik
Protein Families:	Transmembrane

Product images:

Coomassie blue staining of purified TMEM41A protein (Cat# [TP303513]). The protein was produced from HEK293T cells transfected with TMEM41A cDNA clone (Cat# [RC203513]) using MegaTran 2.0 (Cat# [TT210002]).