

Product datasheet for **TP303534M**

ZNF496 (NM_032752) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human zinc finger protein 496 (ZNF496), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203534 protein sequence Red =Cloning site Green =Tags(s)
	MPTALCPRVLAPKESEEPRKMRSPPGENPSPQGELPSPESSRRLFRFRYQEAAGPREALQRLWDLCCGGW LRPERHTKEQILELLVLEQFLAILPREIQSWVRAQEPESGEQAVAAVEALEREPGRPWQWLKHCEDPVVI DDGDSPLDQEQQLPVEPHSDLAKNQDAQPITLAQCLGLPSRPPSQLSGDPVLQDAFLLQEENVRDTQQ V TTLQLPPSRVSPFKDMILCFSEEDWSLLDPAQTGFYGEFIIGEDYGVSMPNDLAAQPDLSQGEENEPRV PELQDLQGKEVPQVSYLDSPLQPFQVEERRKREELQVPEFQACPQTVVPQNTYPAGGNPRSLNSLDEE VTIEIVLSSSGDEDSQHGPYCTEELGSPTEKQRSPLASHRSSTEAGGEVQTSKKSYPVPCNCGKIFRWRVN FIRHLRSRREQEPHECSVCGELFSDSEDLDGHLESHEAQPYRCGACGKSFRLNSHLLSHRRIHLQPD LQPVKEKREQAASEDADKGPKPELNGKAKLSFQCCECGKAFQRHDHLARHRSHFHLKDKARPFQCRYCV K SFTQNYDLLRHERLHMKRRSKQALNSY SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	66.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_116141</u>
Locus ID:	84838
UniProt ID:	<u>Q96IT1</u>
RefSeq Size:	2453
Cytogenetics:	1q44
RefSeq ORF:	1761
Synonyms:	NIZP1; ZFP496; ZKSCAN17; ZSCAN49
Summary:	DNA-binding transcription factor that can both act as an activator and a repressor. [UniProtKB/Swiss-Prot Function]
Protein Families:	Transcription Factors

Product images:



Coomassie blue staining of purified ZNF496 protein (Cat# [TP303534]). The protein was produced from HEK293T cells transfected with ZNF496 cDNA clone (Cat# [RC203534]) using MegaTran 2.0 (Cat# [TT210002]).