

Product datasheet for TP303537

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

SPINT1 (NM_001032367) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human serine peptidase inhibitor, Kunitz type 1 (SPINT1), transcript

variant 3, 20 µg

Species: Human
Expression Host: HEK293T

Expression cDNA Clone >RC203537 protein sequence or AA Sequence: Red=Cloning site Green=Tags(s)

MAPARTMARARLAPAGIPAVALWLLCTLGLQGTQAGPPPAPPGLPAGADCLNSFTAGVPGFVLDTNASVS NGATFLESPTVRRGWDCVRACCTTQNCNLALVELQPDRGEDAIAACFLINCLYEQNFVCKFAPREGFINY LTREVYRSYRQLRTQGFGGSGIPKAWAGIDLKVQPQEPLVLKDVENTDWRLLRGDTDVRVERKDPNQVEL WGLKEGTYLFQLTVTSSDHPEDTANVTVTVLSTKQTEDYCLASNKVGRCRGSFPRWYYDPTEQICKSFVY GGCLGNKNNYLREEECILACRGVQGPSMERRHPVCSGTCQPTQFRCSNGCCIDSFLECDDTPNCPDASDE AACEKYTSGFDELQRIHFPSDKGHCVDLPDTGLCKESIPRWYYNPFSEHCARFTYGGCYGNKNNFEEEQQ CLESCRGISKKDVFGLRREIPIPSTGSVEMAVAVFLVICIVVVVAILGYCFFKNQRKDFHGHHHHPPPTP

ASSTVSTTEDTEHLVYNHTTRPL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW:

Concentration: >0.05 µg/µL as determined by microplate BCA method

53.3 kDa

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by

conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.



SPINT1 (NM_001032367) Human Recombinant Protein - TP303537

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 001027539

Locus ID: 6692 **UniProt ID:** 043278

RefSeq Size: 2342

Cytogenetics: 15q15.1 RefSeq ORF: 1539


Synonyms: HAI; HAI1; MANSC2

Summary: The protein encoded by this gene is a member of the Kunitz family of serine protease

inhibitors. The protein is a potent inhibitor specific for HGF activator and is thought to be involved in the regulation of the proteolytic activation of HGF in injured tissues. Alternative splicing results in multiple variants encoding different isoforms. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Secreted Protein, Transmembrane

Product images:



Coomassie blue staining of purified SPINT1 protein (Cat# TP303537). The protein was produced from HEK293T cells transfected with SPINT1 cDNA clone (Cat# [RC203537]) using MegaTran 2.0 (Cat# [TT210002]).