

Product datasheet for **TP303569L**

HNRPLL (HNRNPLL) (NM_138394) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human heterogeneous nuclear ribonucleoprotein L-like (HNRPLL), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203569 protein sequence Red =Cloning site Green =Tags(s)
	<p>MSSSSSSPRETYEEDREYESQAKRLKTEEGEIDYSAEEGENRREATPRGGGDGGGGGRSFSQPEAGGSHH KVSVPVHVHRLCESVVEADLVEALEKFGTICYVMMMPFKRQALVEFENIDSAKECVTFAADEPVYIAG QQAFFNYSTSKRITRPGNTDDPSGGNKVLLLSIQNPLYPITVDVLYTVCNPVGKVQRIVIFKRNGIQAMV EFESVLCQAQAKAALNGADIYAGCCTLKIEYARPTRLNVIRNDNDSWDYTKPYLGRRDRGKGRQRQAILG EHPSSFRHDGYGSHGPLLPLPSRYRMGSRDTPELVAYPLPQASSYMHGGNPSGSVVMVSGLHQLKMN CS RVFNLFCLYGNIEKVFKMTIPGTALVEMGDEYAVRAVTHLNNVKLFGKRLNVCVSKQHSVPSQIFEL EDGTSSYKDFAMSKNNRFTSAGQASKNIIQPPSCVLHYYNVPLCVTEETFTKLCNDHEVLTFIKYKVFDA KPSAKTSLGLEWECKTDAVEALTALNHYQIRVPNGSNPYTLKLCFSTSSH</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	59.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_612403</u>
Locus ID:	92906
UniProt ID:	<u>Q8WVV9</u>
RefSeq Size:	3053
Cytogenetics:	2p22.1
RefSeq ORF:	1626
Synonyms:	HNRPLL; SRRF
Summary:	HNRNPLL is a master regulator of activation-induced alternative splicing in T cells. In particular, it alters splicing of CD45 (PTPRC; MIM 151460), a tyrosine phosphatase essential for T-cell development and activation (Oberdoerffer et al., 2008 [PubMed 18669861]). [supplied by OMIM, Aug 2008]

Product images:



Coomassie blue staining of purified HNRNPLL protein (Cat# [TP303569]). The protein was produced from HEK293T cells transfected with HNRNPLL cDNA clone (Cat# [RC203569]) using MegaTran 2.0 (Cat# [TT210002]).