

## Product datasheet for **TP303570L**

### VILIP3 (HPCAL1) (NM\_002149) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hippocalcin-like 1 (HPCAL1), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203570 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MGKQNSKLRPEVLQDLRENTFTDHELQEWYKGFLKDCPTGHLTVDEFKKIYANFFPYGDASKFAEHVFR TFDTNGDGTIDFREIIALSVTSRGKLEQKLKWAFSMYDLDDGNGYISRSEMLEIVQAIYKMVSSVMKMPE DESTPEKRTDKIFRQMDTNNDGKLSLEEFIRGAKSDPSIVRLLQCDPSSASQF</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-Myc/DDK
Predicted MW:	22.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_002140</a></u>
Locus ID:	3241
UniProt ID:	<u><a href="#">P37235</a></u>
RefSeq Size:	1779



[View online »](#)

<b>Cytogenetics:</b>	2p25.1
<b>RefSeq ORF:</b>	579
<b>Synonyms:</b>	BDR1; HLP2; VILIP-3
<b>Summary:</b>	The protein encoded by this gene is a member of neuron-specific calcium-binding proteins family found in the retina and brain. It is highly similar to human hippocalcin protein and nearly identical to the rat and mouse hippocalcin like-1 proteins. It may be involved in the calcium-dependent regulation of rhodopsin phosphorylation and may be of relevance for neuronal signalling in the central nervous system. Several alternatively spliced transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Apr 2012]

### Product images:



Coomassie blue staining of purified HPCAL1 protein (Cat# [TP303570]). The protein was produced from HEK293T cells transfected with HPCAL1 cDNA clone (Cat# [RC203570]) using MegaTran 2.0 (Cat# [TT210002]).