

## Product datasheet for TP303572

### C6orf136 (NM\_145029) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 6 open reading frame 136 (C6orf136), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203572 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MEEHLSVMYERLRQELPKLFLQSHDYSLSLDVEFINEILNIRTKGRTWYILSLTLCRFLAWNYFAHLRL EVLQLTRHPENWTLQARWRLVGLPVHLLFLRFYKRDKDEHYRTYDAYSTFYLNSSGLICRHLRDLKLMPSH SPPTPVKLLV GALVALGLSEPEPDLNLCSKP  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	21.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_659466</a>
Locus ID:	221545
UniProt ID:	<a href="#">Q5SQH8</a>




[View online »](#)

RefSeq Size: 1039

Cytogenetics: 6p21.33

RefSeq ORF: 516

### Product images:



Coomassie blue staining of purified C6orf136 protein (Cat# TP303572). The protein was produced from HEK293T cells transfected with C6orf136 cDNA clone (Cat# [RC203572]) using MegaTran 2.0 (Cat# [TT210002]).