

## Product datasheet for TP303606L

### C1orf21 (NM\_030806) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 1 open reading frame 21 (C1orf21), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203606 protein sequence <span style="color: red;">Red</span> =Cloning site <span style="color: green;">Green</span> =Tags(s)  MGCASAKHVATVQNEEEAQKGKNYQNGDVFGEYRIKPVEEVKYMKNGAEEEQKIAARNQENLEKSASS N VRLKTNKEVPGLVHQPRANMHISESQQEFFRMLDEKIEKGRDYCSEEDIT  <span style="color: red;">TR</span> <span style="color: green;">TRPLEQKLISEEDLAANDILDYKDDDDKV</span>
Tag:	C-Myc/DDK
Predicted MW:	13.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_110433</a></u>
Locus ID:	81563
UniProt ID:	<u><a href="#">Q9H246</a></u>
RefSeq Size:	10308


[View online »](#)

Cytogenetics: 1q25.3

RefSeq ORF: 363

Synonyms: PIG13

**Product images:**


Coomassie blue staining of purified C1orf21 protein (Cat# [TP303606]). The protein was produced from HEK293T cells transfected with C1orf21 cDNA clone (Cat# [RC203606]) using MegaTran 2.0 (Cat# [TT210002]).