

## Product datasheet for TP303619

### R3HDM4 (NM\_138774) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 19 open reading frame 22 (C19orf22), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203619 representing NM_138774 Red=Cloning site Green=Tags(s)

MVALENPECGPEAAEGTPGGRRLLPLPSCLPALASSQVKRLSASRRKQHFQAVRNSDLVPAKGRKSL  
QRLNTQYLLTLETDGGLPGLEDGDLAPPASPGIFAEACNNATYVEVWVDFMNRSGEEQERVLYLEDE  
GRSKARRRGPRGEDRRREDPAYTPRECFQRISRRLRAVLKRSRIPMETLETWEERLLRFFSVSPQAVYT  
AMLDNSFERLLLHAVCQYMDLISASADLEGKRQMKVSNRHLDLPPGLLLSAYLEQHS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	30.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_620129</a>
Locus ID:	91300
UniProt ID:	<a href="#">Q96D70</a>



[View online »](#)

RefSeq Size: 1810  
Cytogenetics: 19p13.3  
RefSeq ORF: 804  
Synonyms: C19orf22

**Product images:**

Coomassie blue staining of purified R3HDM4 protein (Cat# TP303619). The protein was produced from HEK293T cells transfected with R3HDM4 cDNA clone (Cat# [RC203619]) using MegaTran 2.0 (Cat# [TT210002]).