

Product datasheet for **TP303635**

RRP9 (NM_004704) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ribosomal RNA processing 9, small subunit (SSU) processome component, homolog (yeast) (RRP9), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203635 protein sequence Red =Cloning site Green =Tags(s)
	<p>MSATAAARKRGKPASGAGAGAGAGKRRRKADSAGDRGKSKGGGKMNEISSDSESESLAPRKPEEEEEEE LEETAQEKKRLAKLYLEQLRQEEEEKAEARAFEEEDQVAGRLKEDVLEQRGRLQKLVAKEIQAPASADIR VLRGHLQSITCLVTPDDSAIFSAAKDCSIKWSVESGRKLHVIPRAKKGAEKPPGHSSHVLCMAISSD GKYLASGDRSKLILIWEAQSCQHLYTFTGHRDAVSGLAFRRGTHQLYSTSHDRSVKVWNVAENSYVETLF GHQDAVAALDALSRECCVTAGGRDGTVRVWVKIPEESQLVFYGHQGSIDCIHLINEEHVMVSGADDGVALW GLSKKRPLALQREAHGLRGEPLGEQPFWISSVAALLNTDLVATGSHSSCVRLWQCCEGFRQLDLLCDIPL VGFINSCLKFSSSGDFLVAGVGQEHRLGRWWRIKEARNSVCIPLRRVPVPPAAGS</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	51.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_004695](#)

Locus ID: 9136

UniProt ID: [O43818](#)

RefSeq Size: 1590


Cytogenetics: 3p21.2

RefSeq ORF: 1425

Synonyms: RNU3IP2; U3-55K

Summary: This gene encodes a member of the WD-repeat protein family. The encoded protein is a component of the nucleolar small nuclear ribonucleoprotein particle (snoRNP) and is essential for 18s rRNA processing during ribosome synthesis. It contains seven WD domains required for nucleolar localization and specific interaction with the U3 small nucleolar RNA (U3 snoRNA). [provided by RefSeq, Oct 2012]

Product images:



Coomassie blue staining of purified RRP9 protein (Cat# TP303635). The protein was produced from HEK293T cells transfected with RRP9 cDNA clone (Cat# [RC203635]) using MegaTran 2.0 (Cat# [TT210002]).