

Product datasheet for TP303658

C1orf76 (FAM163A) (NM_173509) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human family with sequence similarity 163, member A (FAM163A), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203658 protein sequence Red =Cloning site Green =Tags(s)
	 MTAGTWWITGGILATVILLCCIIVLCYCRLQYYCCKKSGTEVADEEEEREHDLPTHPRGPTCNACSSQAL DGRGSLAPLTSEPCSQPCGVAASHCTTCSPPYSSPFYIRTADMVPNGGGGERLSFAPTYYYKEGGPPSLKLA APQSYPTWPGSGREAFTNPRAISTDV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	17.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_775780</u>
Locus ID:	148753
UniProt ID:	<u>Q96GL9</u>



[View online »](#)

RefSeq Size:	2927
Cytogenetics:	1q25.2
RefSeq ORF:	501
Synonyms:	C1orf76; NDSP
Protein Families:	Transmembrane

Product images:

Coomassie blue staining of purified FAM163A protein (Cat# TP303658). The protein was produced from HEK293T cells transfected with FAM163A cDNA clone (Cat# [RC203658]) using MegaTran 2.0 (Cat# [TT210002]).