

## Product datasheet for **TP303665M**

### Seryl tRNA synthetase (SARS) (NM\_006513) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human seryl-tRNA synthetase (SARS), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203665 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MVLDLDFRVDKGGDPALIRETQEKRFPGLVDQLVKADSEWRRRCFRADNLNKLKNLCSKTIGEKMKK  
KEPVGDDSEVPENVLSFDDLTADALANLKVSIKKVRLIDEAILKCD AERIKLEAERFENLREIGNLLH  
PSVPISNDEEDVDNKVERIWGDCTVRKKYSHVDLWMVDGFEGEKGA VAGSRGYFLKGVLFLEQALIQY  
ALRTLGSRGYIPIYTPFFMRKEVMQEVAQLSQFDEELYK VIGKGSEKSDDNSYDEKYLIATSEQPIAALH  
RDEWLRPEDLPIKYAGLSTCFRQEVGSHGRDTRGIFRVHQFEKIEQFVYSSPHDNKSWEMFEEMITTAEE  
FYQSLGIPYHIVNIVSGSLNHAASKKLDLEAWFPGSGAFREL VSCSNCTDYQARRLRIRYGQTKKMMDKV  
EFVHMLNATMCATTCTICAILENYQTEKGITVPEKLKEFMPPGLQELIPFVKPAIEQEPSKKQKQHEG  
SKKKAARDVTLENRLQNMEVTD A

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	58.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online >](#)

RefSeq: [NP\\_006504](#)

Locus ID: 6301

UniProt ID: [P49591](#), [Q0VGA5](#)

RefSeq Size: 1956

Cytogenetics: 1p13.3


RefSeq ORF: 1542

Synonyms: NEDMAS; SARS; SERRS; SERS

**Summary:** This gene belongs to the class II amino-acyl tRNA family. The encoded enzyme catalyzes the transfer of L-serine to tRNA (Ser) and is related to bacterial and yeast counterparts. Multiple alternatively spliced transcript variants have been described but the biological validity of all variants is unknown. [provided by RefSeq, Jul 2010]

**Protein Pathways:** Aminoacyl-tRNA biosynthesis

### Product images:



Coomassie blue staining of purified SARS protein (Cat# [TP303665]). The protein was produced from HEK293T cells transfected with SARS cDNA clone (Cat# [RC203665]) using MegaTran 2.0 (Cat# [TT210002]).