

Product datasheet for TP303681L

GRAMD2B (NM_023927) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human GRAM domain containing 3 (GRAMD3), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203681 protein sequence Red =Cloning site Green =Tags(s)
	<p>MTELQQDVEDTKPAKVLGKRESKLGSAHSEAENGVEEKKKACRSPTAQSPTPSVEADSPDQKKIISLWSK SSFDGASLASDKNDCKTESKNDPKTERKKSSSSSQYKANMHFKLFLSVPTTEPLKQSFTCALQKEILYQ GKLFVSENWICFHSKVFGKDTKISIPAFSVTLIKKTKTALLVPNALIIATVTDRIYFVSLLSRDSTYKLL KSVCGHLENTSVGNPNPSSAENSFRADRPSSLPLDFNDEFSDLDGVVQRRQDMEGYSSSGSQTPESEN SRDFHATESQTVLNVSKGEAKPTRADAHVNRVPEGKAKSLPVQGLSETVGILHKVKSQKCPMLHHILIFY AIVCALIISTFYMYRINTLEEQLGLLTSIVDTHNTEQAAPSGLRSQVQFNVEVLCQELTANIVKLEKI QNNLQKLENGD</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	47.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_076416</u>



[View online »](#)

Locus ID: 65983
UniProt ID: [Q96HH9](#)
RefSeq Size: 2934
Cytogenetics: 5q23.2
RefSeq ORF: 1296
Synonyms: GRAMD3; NS3TP2
Protein Families: Transmembrane

Product images:



Coomassie blue staining of purified GRAMD2B protein (Cat# [TP303681]). The protein was produced from HEK293T cells transfected with GRAMD2B cDNA clone (Cat# [RC203681]) using MegaTran 2.0 (Cat# [TT210002]).