

Product datasheet for TP303687M

TADA3L (TADA3) (NM_006354) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human transcriptional adaptor 3 (NGG1 homolog, yeast)-like (TADA3L), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203687 protein sequence Red =Cloning site Green =Tags(s)
	MSELKDCPLQFHDFKSV DHLKVCPRYTAVLARSEDDGIGIEELDTLQLELETLSSASRRLRVLEAETQI LTDWQDKKGD RRFLKLRDHELGAPPKHGKPKKQKLEGKAGHGPGPGPRPKSKNLQPKIQEYFTDDPI DVPRIPKNDAPNRFWASVEPYCADITSEEVRTLEELLKPPEDA EAHYKIPPLGKHYSQRWAQEDLLEEQK DGARAAAVADKKKGLMGPLTELDTKDVDALLKSEAQHEQPEDGCPFGALTQ RLLQALVEENIISPMEDS PIPDMSGKESGADGASTSPRNQNKPFVPHTKSLESRIKEELIAQGLLESEDRPAEDSEDEVLAE LRKRQ AELKALSAHNRTKKHDLRLAKEEVS RQELRQRV RMADNEVMDAFRKIMAA RQKKRTPTKKEKDQAWKTL KERESILKLLDG
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	48.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_006345](#)

Locus ID: 10474

UniProt ID: [O75528](#), [A8K899](#)

RefSeq Size: 2530

Cytogenetics: 3p25.3


RefSeq ORF: 1296

Synonyms: ADA3; hADA3; NGG1; STAF54; TADA3L

Summary: DNA-binding transcriptional activator proteins increase the rate of transcription by interacting with the transcriptional machinery bound to the basal promoter in conjunction with adaptor proteins, possibly by acetylation and destabilization of nucleosomes. The protein encoded by this gene is a transcriptional activator adaptor and a component of the histone acetyl transferase (HAT) coactivator complex which plays a crucial role in chromatin modulation and cell cycle progression. Along with the other components of the complex, this protein links transcriptional activators bound to specific promoters, to histone acetylation and the transcriptional machinery. The protein is also involved in the stabilization and activation of the p53 tumor suppressor protein that plays a role in the cellular response to DNA damage. Alternate splicing results in multiple transcript variants of this gene. [provided by RefSeq, May 2013]

Protein Families: Transcription Factors

Product images:



Coomassie blue staining of purified TADA3 protein (Cat# [TP303687]). The protein was produced from HEK293T cells transfected with TADA3 cDNA clone (Cat# [RC203687]) using MegaTran 2.0 (Cat# [TT210002]).