

Product datasheet for **TP303724**

CENTB1 (ACAP1) (NM_014716) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ArfGAP with coiled-coil, ankyrin repeat and PH domains 1 (ACAP1), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203724 protein sequence Red =Cloning site Green =Tags(s)
	<p>MTVKLDFEELKDSRFRASIELVEAEVSELETRLEKLLKLGTLLESGRHYLAASRAFWGICDLARLG PPEPMMAECLKFTVSLNHKLDSHAELLDATQHTLQQQIQTLVKEGLRGFREARRDFWRGAESLEAALTH NAEVRPRRAQEAEAGAAALRTARAGYRGRALDYALQINVIEDKRKFDIMEFVLRRLVEAQATHFQQGHEEL SRLSQYRKELGAQLHQLVLNSAREKRDMEQRHVLLKQKELGGEEPEPSLREGPGGLVMEGHLFKRASNAF KTWSRRWFITQSNQLVYQKKYKDPVTVVWDDLRLCTVKLCPDSERRFCFEWVSTSKSCLLQADSERLLQL WWSAVQSSIASAFSQARLDDSPRGGQSGHLAGSAATLGS GGMARGREPGGVGHVVAQVQSV DGNAC CDCREPAPEWASINLGVTLICQCSGIHRSLGVHFSKVRSLTLDSEPELVKLMCELGNVIINQIYEARVE AMAVKKPGPSCSRQEKEAWIHAKYVEKKFLTKLPEIRRRGGRRGRPRGQPPVPPKPSIRPRPGSLRSKPE PPSEDLGSLHPGALLFRASGHPPSLPTMADALAHGADVNVWNGGQDNATPLIQATAANSL LACEFLLQNG ANVNQADSAGRGPLHHATILGHTGLACLFLKRGADLGARDSEGRDPLTIAMETANADIVTLLRLAKMREA EAAQGQAGDETYLDIFRDFSLMASDDPEKLSRRSHDLHTL</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	81.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_055531
Locus ID:	9744
UniProt ID:	Q15027
RefSeq Size:	2523
Cytogenetics:	17p13.1
RefSeq ORF:	2220
Synonyms:	CENTB1
Summary:	GTPase-activating protein (GAP) for ADP ribosylation factor 6 (ARF6) required for clathrin-dependent export of proteins from recycling endosomes to trans-Golgi network and cell surface. Required for regulated export of ITGB1 from recycling endosomes to the cell surface and ITGB1-dependent cell migration.[UniProtKB/Swiss-Prot Function]
Protein Pathways:	Endocytosis

Product images:



Coomassie blue staining of purified ACAP1 protein (Cat# TP303724). The protein was produced from HEK293T cells transfected with ACAP1 cDNA clone (Cat# [RC203724]) using MegaTran 2.0 (Cat# [TT210002]).