

Product datasheet for TP303737L

UCHL5IP (HAUS7) (NM_017518) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human UCHL5 interacting protein (UCHL5IP), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203737 protein sequence Red =Cloning site Green =Tags(s)

MGGARLGARNMAGQDAGCGRGGDDYSEDEGDSSVSRAAVEVFGKLDLNCPLFLEGLYITEPKTIQELLCS
PSEYRLEILEWMCTRVWPSLQDRFSSLKGVPTVEVKIQEMTKLGHMLCAPDDQELLKGCACAQKQLHFM
DQLLDTIRSLTIGSSCSSLMEHFEDTREKNEALLGELFSSPHLQMLLNPECDPWPLDMQPLLNKQSDDW
QWASASAKSEEEEEKLAELARQLQESA AKLHALRTEYFAQHEQGAAAGAADISTLDQKRLRLVTSDFHQLIL
AFLQVYDDELGECCQRPGLHPCGPPIIQATHQNLTSYSQLLQVVMVAADTSAKAVETVKKQOGEQICWG
GSSVMSLATKMNELMEK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	40.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_059988</u>




[View online »](#)

Locus ID: 55559
UniProt ID: [Q99871](#)
RefSeq Size: 1378
Cytogenetics: Xq28
RefSeq ORF: 1104
Synonyms: UCL5IP; UIP1

Summary: This gene encodes a subunit of the augmin complex, which regulates centrosome and mitotic spindle integrity, and is necessary for the completion of cytokinesis. The encoded protein was identified by interaction with ubiquitin C-terminal hydrolase 37. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2012]

Product images:



Coomassie blue staining of purified HAUS7 protein (Cat# [TP303737]). The protein was produced from HEK293T cells transfected with HAUS7 cDNA clone (Cat# [RC203737]) using MegaTran 2.0 (Cat# [TT210002]).