

## Product datasheet for **TP303800**

### EFS (NM\_032459) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human embryonal Fyn-associated substrate (EFS), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203800 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MAIATSVYVPPPARPCPTSGPPAGPCPPSPDLIYKIPRASGTQLAAPRDALEVYDVPPTALRVPSSGPY DCPASFSHPLTRVAPQPPGEDDAPYDVPLTPKPPAELEPDLEWEGGREPGPIYAAPSNLKRASALLNLY EAPPELLADGEGGGTDEGIYDVPLLGPEAPPSPEPPGALASHDQDTLAQLLARSPPPPHRPRLPESAESLS RRPLPALPVPEAPSPSPVPSPAPGRKGSIQDRPLPPPPRLPGYGGPKVEGDPEGREMEDDPAGHHNEYE GIPMAEEYDYVHLKGMDKAQGSRPPDQACTGDPELPERGMPAPQEALSPGEPLVSTGDLQLLYFYAGQ C QSHYSALQAAVAALMSSTQANQPPRLFVPHSKRVVAAHRLVFVGDTLGRLAASAPLRAQVRAAGTALG Q ALRATVLAVKGAALGYPSIPAQEMVQCVELAGQALQFTLLTSLAP</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-Myc/DDK
Predicted MW:	48.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_115835</a>
<b>Locus ID:</b>	10278
<b>UniProt ID:</b>	<a href="#">O43281</a>
<b>RefSeq Size:</b>	2851
<b>Cytogenetics:</b>	14q11.2
<b>RefSeq ORF:</b>	1404
<b>Synonyms:</b>	CAS3; CASS3; EFS1; EFS2; HEFS; SIN
<b>Summary:</b>	The protein encoded by this gene is a member of the CAS (CRK-associated substrate) family of adaptor proteins which typically serve as scaffolds for the assembly of larger signaling complexes. These complexes form at the cell surface where integrin binding leads to the subsequent phosphorylation of a CAS protein. Additional binding of SRC family kinases leads to CAS hyperphosphorylation and the creation of binding sites for CRK and other proteins that cause actin cytoskeleton reorganization. This gene plays a role in integrin-mediated cell attachment, spreading, and migration and also plays a role in both normal and malignant cellular transformation. This broadly expressed gene has been shown to play a role in neurite outgrowth and its expression in the thymus and lymphocytes is important for T cell maturation and the development of immunological self-tolerance. Alternative splicing of this gene results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2020]

## Product images:



Coomassie blue staining of purified EFS protein (Cat# TP303800). The protein was produced from HEK293T cells transfected with EFS cDNA clone (Cat# [RC203800]) using MegaTran 2.0 (Cat# [TT210002]).