

Product datasheet for TP303818L

Integrin beta 1 (ITGB1) (NM_002211) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human integrin, beta 1 (fibronectin receptor, beta polypeptide, antigen CD29 includes MDF2, MSK12) (ITGB1), transcript variant 1A, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203818 protein sequence Red=Cloning site Green=Tags(s)

MNLQPIFWIGLISSVCCVFAQTDENRCLKANAKSCGECIQAGPNCGWCTNSTFLQEGMPTSARCDDLEAL
 KKKGCPPDDIENPRGSKDIKKNKNVTNRSKGTAEKLPEDITQIQPQQLVLRRLRSGEPQTFTLKFRAED
 YPIDLYYLMDSL SYSMKDDLENVKS LGTDL MNEMRRITSDFRIGFGSFVEKTVMPYISTTPAKLRNPCTSE
 QNCTSPFSYKNVLSLTNKGEVFNELVGKQRISGNLDSPEGGFDAIMQVAVCGSLIGWRNVTRLLVFSTDA
 GFHFAGDGKLGIVLPNDGQCHLENNMYTMSHYDYPSIAHLVQKLSENNIQTIFAVTEEFQPVYKELKN
 LIPKSAVGTLSANSSNVIQLIIDAYNSLSSEVILENGKLESGVTISYKSYCKNGVNGTGENGGRKCSNISI
 GDEVQFEISITSNKCPKSDSDFKIRPLGFTEEVILQYICECECQSEGIPESPKCHEGNGTFECGACR
 CNEGRVGRHCECSTDEVNSEMDAYCRKENSSEICSNNGECVCGQCVCRRKRDNTNEIYSGKFCECDNFNC
 DRSNGLICGGNGVCKCRVCECNPNYTGACDCSLDTSTCEASNGQICNGRGICECGVCKCTDPKFQGGQTC
 EMCQTCLGVCAEHKECVQCRAFNGEKKDTCTQECSYFNITKVESRDKLPQPVQDPVSHCKEKDVDDCW
 FYFTYSVNGNNEVMVHVVENPECPTGPDIIPIVAGVAGIVLIGLALLLIWKLLMIIHDRREFAKFEKEK
 MNAKWDTGENPIYKSAVTTVNPKEYEGK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	86.1 kDa
Concentration:	>0.1 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_002202
Locus ID:	3688
UniProt ID:	P05556
RefSeq Size:	3879
Cytogenetics:	10p11.22
RefSeq ORF:	2394
Synonyms:	CD29; FNRB; GPIIA; MDF2; MSK12; VLA-BETA; VLAB
Summary:	Integrins are heterodimeric proteins made up of alpha and beta subunits. At least 18 alpha and 8 beta subunits have been described in mammals. Integrin family members are membrane receptors involved in cell adhesion and recognition in a variety of processes including embryogenesis, hemostasis, tissue repair, immune response and metastatic diffusion of tumor cells. This gene encodes a beta subunit. Multiple alternatively spliced transcript variants which encode different protein isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane
Protein Pathways:	Arrhythmogenic right ventricular cardiomyopathy (ARVC), Axon guidance, Cell adhesion molecules (CAMs), Dilated cardiomyopathy, ECM-receptor interaction, Focal adhesion, Hypertrophic cardiomyopathy (HCM), Leukocyte transendothelial migration, Pathogenic Escherichia coli infection, Pathways in cancer, Regulation of actin cytoskeleton, Small cell lung cancer

Product images:



Coomassie blue staining of purified ITGB1 protein (Cat# [TP303818]). The protein was produced from HEK293T cells transfected with ITGB1 cDNA clone (Cat# [RC203818]) using MegaTran 2.0 (Cat# [TT210002]).