

Product datasheet for **TP303835**

PSMG3 (NM_032302) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human proteasome (prosome, macropain) assembly chaperone 3 (PSMG3), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203835 protein sequence Red=Cloning site Green=Tags(s)

MEDTPLVISKQKTEVWCGVPTQWVCTAFSSHILVWTQFGKMGTLVSLEPSSVASDVSKPVLTTKVLGQ
DEPLIHVFAKNLVAFVSQEAGNRAVLLAVAVKDKSMEGLKALREVIRVCQVW

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	12.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_115678</u>
Locus ID:	84262
UniProt ID:	<u>Q9BT73</u> , <u>A0A024R806</u>
RefSeq Size:	1395



[View online »](#)

Cytogenetics:	7p22.3
RefSeq ORF:	366
Synonyms:	C7orf48; PAC3
Summary:	Chaperone protein which promotes assembly of the 20S proteasome. May cooperate with PSMG1-PSMG2 heterodimers to orchestrate the correct assembly of proteasomes. [UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified PSMG3 protein (Cat# TP303835). The protein was produced from HEK293T cells transfected with PSMG3 cDNA clone (Cat# [RC203835]) using MegaTran 2.0 (Cat# [TT210002]).