

Product datasheet for **TP303848**

TIM 4 (TIMD4) (NM_138379) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human T-cell immunoglobulin and mucin domain containing 4 (TIMD4), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203848 protein sequence Red =Cloning site Green =Tags(s)

MSKEPLILWLMIEFWWLYLTPVTSETVTEVLGHRVTLPCLYSSWSHNSNSMCWKGKQCPYSGCKEALIR
TDGMRVTSRKSAYRLQGTIPRGDVSLTILNPSESDSGVYCCRIEVPGWFNVDVKINVRLNLQRASSTHTR
TATTTTRTTTTPTTTRQMTTTPAALPTTWVTPDLTTGTPLQMTTIAVFTTANTCLSLTPSTLPEEAT
GLLTPEPSKEGPILTAESVTLPSDSWSSAESTSADTVLLTSKESKVDLPSTSHVSMWKTSDSVSSPQP
GASDTAVPEQNKTTKTGQMDGIPMSMKNEMPISQLLMIIAPSLGFVLFALFVAFLLRGKLMETYCSQKHT
RLDYIGDSKNVLDVQHGREDDEDGLFTL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	41.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_612388



[View online »](#)

Locus ID:	91937
UniProt ID:	Q96H15
RefSeq Size:	1374
Cytogenetics:	5q33.3
RefSeq ORF:	1134
Synonyms:	SMUCKLER; TIM4
Summary:	Phosphatidylserine receptor that enhances the engulfment of apoptotic cells. Involved in regulating T-cell proliferation and lymphotoxin signaling. Ligand for HAVCR1/TIMD1 (By similarity).[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome, Transmembrane

Product images:



Coomassie blue staining of purified TIMD4 protein (Cat# TP303848). The protein was produced from HEK293T cells transfected with TIMD4 cDNA clone (Cat# [RC203848]) using MegaTran 2.0 (Cat# [TT210002]).