

Product datasheet for **TP303873M**

Rab5 (RAB5A) (NM_004162) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human RAB5A, member RAS oncogene family (RAB5A), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203873 protein sequence Red =Cloning site Green =Tags(s)
	<p>MASRGATRPNGPNTGNKICQFKLVLLGESAVGKSSLVLRVFKGQFHEFQESTIGAAFLTQTVCLDDTTVK FEIWDTAGQERYHSLAPMYRGAQAAIVVYDITNEESFARAKNWWKELQRQASPNIVIALSGNKADLANK RAVDFQEAQSYADDNSLLFMETSAKTSMNVNEIFMAIAKKLPKNPQNPGANSARGRGVDLTEPTQPTRN QCCSN</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	23.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_004153
Locus ID:	5868
UniProt ID:	P20339 , A0A024R2K1



[View online »](#)

RefSeq Size: 2548

Cytogenetics: 3p24.3

RefSeq ORF: 645


Synonyms: RAB5

Summary: The small GTPases Rab are key regulators of intracellular membrane trafficking, from the formation of transport vesicles to their fusion with membranes. Rabs cycle between an inactive GDP-bound form and an active GTP-bound form that is able to recruit to membranes different sets of downstream effectors directly responsible for vesicle formation, movement, tethering and fusion. RAB5A is required for the fusion of plasma membranes and early endosomes (PubMed:10818110, PubMed:14617813, PubMed:16410077, PubMed:15378032). Contributes to the regulation of filopodia extension (PubMed:14978216). Required for the exosomal release of SDCBP, CD63, PDCD6IP and syndecan (PubMed:22660413). Regulates maturation of apoptotic cell-containing phagosomes, probably downstream of DYN2 and PIK3C3 (By similarity).[UniProtKB/Swiss-Prot Function]

Protein Families: Druggable Genome

Protein Pathways: Amyotrophic lateral sclerosis (ALS), Endocytosis

Product images:



Coomassie blue staining of purified RAB5A protein (Cat# [TP303873]). The protein was produced from HEK293T cells transfected with RAB5A cDNA clone (Cat# [RC203873]) using MegaTran 2.0 (Cat# [TT210002]).