

Product datasheet for **TP303894M**

CD147 (BSG) (NM_198589) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human basigin (Ok blood group) (BSG), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203894 protein sequence Red =Cloning site Green =Tags(s)
	MAAALFVLLGFALLGTHGASGAAGTVFTTVEDLGSKILLTCSLNDSEVTGHRWLKGGWVKEDALPGQ KTEFKVDSDDQWGEYSCVFLPEPMGTANIQLHGPPRVKAVKSSEHINEGETAMLVCKSESVPVTDWAWY KITDSEDKALMNGSESRRFFVSSSQGRSELHIENLNMEADPGQYRCNGTSSKGSQDAITLRVRSHLAALW PFLGIVAEVLVLVTIIFIYEKRRKPEDVLDDDDAGSAPLKSSGQHQNCKGKNVRQRNSS
	TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	27.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_940991</u>
Locus ID:	682
UniProt ID:	<u>P35613</u> , <u>Q54A51</u>



[View online »](#)

RefSeq Size:	1759
Cytogenetics:	19p13.3
RefSeq ORF:	807
Synonyms:	5F7; CD147; EMMPRIN; EMPRIN; HAb18G; OK; SLC7A11; TCSF
Summary:	The protein encoded by this gene, basigin, is a plasma membrane protein that is important in spermatogenesis, embryo implantation, neural network formation, and tumor progression. Basigin is also a member of the immunoglobulin superfamily, ubiquitously expressed in various tissues. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2020]
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane

Product images:



Coomassie blue staining of purified BSG protein (Cat# [TP303894]). The protein was produced from HEK293T cells transfected with BSG cDNA clone (Cat# [RC203894]) using MegaTran 2.0 (Cat# [TT210002]).