

Product datasheet for TP303911

AMPK beta 1 (PRKAB1) (NM_006253) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human protein kinase, AMP-activated, beta 1 non-catalytic subunit (PRKAB1), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203911 protein sequence Red =Cloning site Green =Tags(s) MGNTSSERAALERHGGHKTPRRDSSGGTKDGRPKILMDSPEDADLFHSEEIKAPEKEEFLAWQHDLEVN DKAPAQARPTVFRWTGGGKEVYLSGSFNNWSKLPLTRSHNNFVAILDLPEGEHQYKFFVDGQWTHDPSE P IVTSQLGTVNNIIQVKKTDFEVDALMVDSQKCSDVSELSSPPGPYHQEPYVCKPEERFRAPPILPPHL LQVILNKDTGISCDPALLPEPNHVMLNHLYALSIDKGVMVLSATHRYKKKYVTLLYKPI TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	30.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_006244</u>
Locus ID:	5564


[View online »](#)

UniProt ID:	Q9Y478
RefSeq Size:	2412
Cytogenetics:	12q24.23
RefSeq ORF:	810
Synonyms:	AMPK; HAMPKb
Summary:	The protein encoded by this gene is a regulatory subunit of the AMP-activated protein kinase (AMPK). AMPK is a heterotrimer consisting of an alpha catalytic subunit, and non-catalytic beta and gamma subunits. AMPK is an important energy-sensing enzyme that monitors cellular energy status. In response to cellular metabolic stresses, AMPK is activated, and thus phosphorylates and inactivates acetyl-CoA carboxylase (ACC) and beta-hydroxy beta-methylglutaryl-CoA reductase (HMGCR), key enzymes involved in regulating de novo biosynthesis of fatty acid and cholesterol. This subunit may be a positive regulator of AMPK activity. The myristoylation and phosphorylation of this subunit have been shown to affect the enzyme activity and cellular localization of AMPK. This subunit may also serve as an adaptor molecule mediating the association of the AMPK complex. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome
Protein Pathways:	Adipocytokine signaling pathway, Hypertrophic cardiomyopathy (HCM), Insulin signaling pathway

Product images:



Coomassie blue staining of purified PRKAB1 protein (Cat# TP303911). The protein was produced from HEK293T cells transfected with PRKAB1 cDNA clone (Cat# [RC203911]) using MegaTran 2.0 (Cat# [TT210002]).